


COMMONLY USED RETCAM COUPLING MEDIA



GENTEAL GEL (NOVARTIS)
USA, Canada, Australia, Singapore




TEAR-GEL (NOVARTIS)
Canada




VISCOTEAR (NOVARTIS)
United Kingdom
Switzerland



VIDISIC (NOVARTIS)
Germany, Netherlands, Russia, Belarus,
China, Singapore, Romania



POLY GEL (ALCON)
Australia



LIPOSIC (NOVARTIS)
UAE, South Korea



LACRINORM
France, Iran



OTHER COUPLING MEDIA

CORNEREGEL (5MG)
Russia, Belarus, Lithuania

ALCON DRY EYE GEL - OOGGEL (10G)
Netherlands

XYLOCAINE 2% JELLY
India (Chennai)

OFLOXACIN EYE OINTMENT (Antibiotic Gel)
Shanghai, China
Seems to work well as is quite sticky and not a great amount is necessary.

HYDROXYPROPYL METHYLCELLULOSE (Ophthalmic)
Used for general gonioscopy but not always viscous enough to provide adequate coupling for peripheral retinal or iridocorneal imaging using the RetCam.



LACRYVISC (ALCON)
Singapore



THILO-TEARS GEL (ALCON)
Romania



VISTAGEL
Peru